# **GSD NovaType Select L452R SARS-CoV-2 RT-PCR**

**RUO**

**FOR RESEARCH USE ONLY  
Not for use in diagnostic procedures**

|  |    |
|--|----|
| English .....  | 2  |
| Deutsch .....  | 7  |
| Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia .....              | 15 |
| Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas .....          | 15 |
| Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos ..... | 16 |

---

|                 |            |                         |
|-----------------|------------|-------------------------|
| Product Number: | PCOVL452R3 | (1 x 96 Determinations) |
|                 | PCOVL452R4 | (4 x 96 Determinations) |

---

## ENGLISH

### 1. INTRODUCTION

End of 2019, a novel respiratory disease emerged in the city of Wuhan, Hubei Province of the People's Republic of China, and soon spread rapidly within the country and worldwide. The causative agent was identified as severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). SARS-CoV-2 (2019-nCoV), like the closely related SARS coronavirus (SARS-CoV), belongs to the genus Betacoronavirus within the family of coronaviruses. The zoonotic reservoir of the virus appears to be bats.

Coronaviruses are enveloped, positive single-stranded large RNA viruses that infect humans, but also a wide range of animals. The common human coronaviruses NL63, 229E, OC43 and HKU1 are widespread especially throughout the winter months. They are responsible for up to one third of all acute respiratory diseases, typically with mild symptoms (common cold). More than 80 % of the adult population have antibodies against human coronaviruses. The immunity from previous infections lasts only for a short period of time. Therefore, reinfections with the same pathogen are possible just after one year.

SARS-CoV-2 is predominantly transmitted by droplet infection via coughing or sneezing and through close contact with infected patients. In theory, smear infection and infection through the conjunctiva of the eyes are also possible.

The incubation period is in the median 5–6 days (and up to 14 days maximum).

The clinical manifestations of SARS-CoV-2-related COVID-19 disease include fever, cough, respiratory problems and fatigue. In most patients the infection manifests with symptoms of a mild febrile illness with irregular lung infiltrates.

The initial clinical sign of COVID-19 which allowed case detection was pneumonia. But it turned out that the course of the disease is non-specific and varies widely, from asymptomatic courses to severe pneumonia with lung failure and death. However, based on current knowledge, around 80 % of the illnesses are mild to moderate.

Although severe courses of the disease also occur in younger patients and people without previous illness, the following groups of people have an increased risk of serious forms of the disease: elderly people (with a steadily increasing risk from around 50-60 years of age), smokers and people with certain diseases of the cardiovascular system or the lungs, patients with chronic liver diseases, diabetes mellitus, cancer, or patients with a weakened immune system (e.g. due to immune deficiencies or by taking drugs that suppress the immune system).

Treatment focuses on supportive measures depending on the severity of the clinical picture (e.g. oxygen administration, fluid balance management, etc.) as well as the treatment of relevant underlying diseases. Various specific therapeutic approaches (directly antiviral effective, immunomodulatory effective) have been and are being investigated in studies during the course of the SARS-CoV-2 pandemic.

The SARS-CoV-2 trimeric spike glycoprotein (S) is exposed on the surface of virions and mediates viral entry into host cells. The S protein constitutes the primary target of all current leading vaccine candidates. This large, exposed protein is readily targeted by neutralizing antibodies, which indirectly creates selective pressure for the emergence of evasion mutations. The propensity for S to mutate may limit its future use in serological assays and vaccines, as antibodies directed against the current variant may not bind emerging mutated epitopes.

Gene mutations in emerging SARS-CoV-2 variants can affect transmission, virulence, and immunity. Variants for which there is evidence (variants of concern, VOC) or suspicion (variants of interest, VOI) of increased transmissibility, more severe disease (e.g., increased hospitalizations or deaths), significant reduction in neutralization by antibodies, reduced efficacy of treatments or vaccines, or diagnostic detection errors are considered to be under special surveillance worldwide. Current VOCs and VOIs include the lineages, B.1.1.7 (501Y.V1), B.1.351 (20H/501Y.V2), B.1.1.28 (P.1; 20J/501Y.V3), B.1.617 (G/452.V3) and B.1.427/B.1.429 (CAL.20C). They all exhibit mutations in the spike protein that are associated with increased transmissibility and/or suspected to negatively affect vaccine efficacy by avoiding the immune response. The lineages B.1.617 (G/452.V3) and B.1.427/B.1.429 (CAL.20C) are both found to be carrier of the mutation L452R in addition to other mutations. This mutation is suspected to be associated with both increased infectivity and reduced immune response and is therefore of relevance to the global pandemic situation.

| Species   | Disease  | Symptoms e.g.  | Transmission route   |
|---|----------|--|--|
| SARS-CoV-2<br>(severe acute respiratory syndrome coronavirus 2) | COVID-19 | the course of the disease is unspecific, diverse and varies greatly, from asymptomatic courses to severe pneumonia with lung failure and death | primary mode of transmission: droplet infection;<br>smear infections and infections via the conjunctiva of the eyes are theoretically possible |

Infection or presence of pathogen may be identified by

- Nucleic acid Amplification Techniques (NAT): e.g. RT-PCR
- Serology: detection of antibodies by e.g. ELISA

### 2. INTENDED USE (FOR RESEARCH USE ONLY!)

The GSD NovaType Select L452R SARS-CoV-2 is a research use only (RUO) real-time reverse transcription PCR (RT-PCR) assay designed for the qualitative detection of the SARS-CoV-2 (Severe acute respiratory syndrome coronavirus 2) S gene mutation L452R. Sample material should be genomic RNA extracted from human upper respiratory (nasal wash/smear, nasopharyngeal wash/smear, oropharyngeal swab) specimen types that have previously tested positive for SARS-CoV-2 RNA by diagnostic RT-PCR methods e.g. GSD NovaPrime® SARS-CoV-2 (COVID-19). The GSD NovaType Select L452R SARS-CoV-2 is expressly not intended for diagnosis of SARS-CoV-2 infection and not for use in diagnostic procedures. It is only designed to aid investigations related to the prevalence and frequency of virus variants.

### 3. PRINCIPLE OF THE ASSAY

The qualitative determination of specific RNA is based on Real-Time reverse-transcription Polymerase Chain Reaction (RT-PCR) technology. The kit contains specific primers and probes labelled with fluorescent reporter and quencher molecules for amplification and simultaneous detection of the specific RNA sequence representing the **SARS-CoV-2 S gene mutation L452R**. The S gene specific probe is labelled with the fluorophore FAM™. As internal process control, the human Ribonuclease P (RNase P) is detected using an CY5™ labelled probe. The human RNase P gene primer and probe set serves as an internal positive control to monitor sample quality, RNA extraction, and for detection of inhibitors of the PCR reaction.

The GSD NovaType Select L452R SARS-CoV-2 has been validated for following real-time PCR platforms:

- AriaMx™ (Agilent Technologies)
- AriaDx™ (Agilent Technologies)

### 4. MATERIALS

#### 4.1. Reagents Provided

| Cap         | Symbol          | Component   | PCOVL452R3 (PCOVL452R4) |                 |
|-------------|-----------------|---|-------------------------|-----------------|
|             |                 |   | Volume per Vial [µL]    | Number of Vials |
| green       | <b>E-MIX</b>    | RT-PCR Enzyme Mix<br>(reverse transcriptase, hot-start DNA polymerase, RNase inhibitor, nucleotides, magnesium, stabilizers and buffer) | 550                     | 1 (4)           |
| blue/white  | <b>PP L452R</b> | Primer-Probe-Mix<br>(for detection of mutation L452R and for detection of RNase P gene)   | 350                     | 1 (4)           |
| red         | <b>PC</b>       | Positive Control<br>(plasmid DNA representing SARS-CoV-2 S gene mutation L452R)   | 150                     | 1 (4)           |
| transparent | <b>NFW</b>      | Nuclease Free Water   | 500                     | 1 (4)           |

#### 4.2. Materials and Equipment needed, but not provided

- Deionized nuclease-free and nucleic-acid free water
- Real-Time PCR instrument, equipped with suitable set of optical cartridges (for already validated instruments refer to 3. PRINCIPLE OF THE ASSAY; for dyes used refer to 7.2.)  
Alternative Real-Time PCR instruments might also be appropriate. Their suitability for use with the GSD NovaType Select L452R SARS-CoV-2 has to be validated by the user.
- Appropriate Real-Time PCR consumables (e.g. disposable tubes, reaction plates, corresponding optical closing materials)
- Benchtop microcentrifuge
- Centrifuge with a rotor for microtiter plates
- Vortex mixer
- Adjustable pipettes in relation to reaction setup
- Disposable DNase/RNase free pipette tips with filters
- Disposable powder-free gloves

### 5. STABILITY AND STORAGE

The GSD NovaType Select L452R SARS-CoV-2 kit is shipped on dry ice and all components should arrive frozen.

- All components have to be stored between -30 °C and -15 °C immediately after arrival.
- Repeated freeze thaw cycles (more than two) of reagents should be avoided, since this might affect the performance of the kit. Reagents should be frozen in aliquots if they are used intermittently.
- Keep unfrozen storage (e.g. storage on ice) as short as possible.
- Protect **E-MIX** and **PP L452R** from light and keep frozen until use.

### 6. SAMPLE PREPARATION

- Extracted RNA or total nucleic acid extracted from human respiratory specimen types (nasal lavage/smear, nasopharyngeal lavage/smear, oropharyngeal swab, bronchoalveolar lavage) previously tested positive for SARS-CoV-2 RNA by RT-PCR methods e.g. GSD NovaPrime® SARS-CoV-2 (COVID-19) is the starting material for the GSD NovaType Select L452R SARS-CoV-2.
- The quality and the quantity of the extracted RNA has a crucial effect on the performance of the entire RT-PCR test system. Make sure that the nucleic acid extraction method is compatible with Real-Time PCR technology.

## 7. ASSAY PROCEDURE

### 7.1. Reaction Setup

- Please read the instructions for use carefully **before** performing the assay. Reliability of results depends on following strictly the instructions for use.
- Before use make sure that all samples and reagents are thawed completely, mixed by up and down pipetting or vortexing and centrifuged briefly.
- **Important:** **E-MIX** is viscous. Briefly spin down the tube to ensure material has not lodged in the cap or side of tube. Be sure to pipette and dispense carefully, and use pipette tips suitable for pipetting viscous liquids.
- The use of **NFW** as no template control (NTC) is highly recommended.
- Define the positions (wells) for samples and controls (**PC** or **NFW**) on the plate.

| Reaction Setup                    |              |
|-----------------------------------|--------------|
| Component                         | Volume       |
| <b>E-MIX</b>                      | 5 µL         |
| <b>PP L452R</b>                   | 3 µL         |
| Sample or <b>PC</b> or <b>NFW</b> | 12 µL        |
| <b>Total Volume</b>               | <b>20 µL</b> |

- Close the optical reaction plate with corresponding optical closing material.
- Centrifuge the optical reaction plate in a centrifuge with a rotor for microtiter plates for 60 seconds at approximately 1000 x g (~ 3000 rpm).

### 7.2. Programming the Real-Time PCR Instrument

Regarding setup and programming of the Real-Time PCR instrument, please use the manual of the respective instrument.

For detailed programming instructions regarding the use of the GSD NovaType Select L452R SARS-CoV-2 on specific Real-Time PCR instruments please contact NovaTec Immundiagnostica GmbH.

| RT-PCR Run Settings |       |
|---------------------|-------|
| Reaction Volume     | 20 µL |
| Passive Reference   | None  |

#### Fluorescent Detectors / Fluorophores

Detection of the amplified viral nucleic acid fragment for the S gene mutation L452R is performed in the FAM™ detection channel with MGBEQ Quencher. Detection of the amplified internal process control RNase P is performed in the detection channel CY5™ with BHQ2 Quencher.

| Target       | Fluorophore/Detection Channel | Quencher |
|--------------|-------------------------------|----------|
| S gene L452R | FAM™                          | MGBEQ    |
| RNase P      | CY5™                          | BHQ2     |

#### Temperature Profile and Data Collection

| No. of Cycles | Temperature | Time (min) | Data Collection                                    |
|---------------|-------------|------------|--|
| 1             | 50 °C       | 10:00      | -  |
| 1             | 95 °C       | 03:00      | -  |
| 40            | 95 °C       | 00:10      | -  |
|               | 60 °C       | 00:30      | Fluorescence measurement at the end of every cycle |

Before starting the test run, please check the settings for cycles, temperature and time.

## 8. RESULTS

Data analysis should be performed with the software of the used real-time PCR device according to manufacturer's instructions.

Analysis settings:

| Setting   | Recommendation   |
|-----------|--|
| Threshold | Enter a value for the threshold so that the threshold is: <ul style="list-style-type: none"> <li>Above the background</li> <li>Below the plateau and linear regions of the amplification curve</li> <li>Within the exponential phase of the amplification curve</li> </ul> |
| Baseline  | Select start and end cycle values so that possible initial signal noise is not included and the baseline ends before significant fluorescence is detected.   |

### 8.1. Interpretation of Results

The test run is valid only if the RT-PCR run is complete.

The GSD NovaType Select L452R SARS-CoV-2 test protocol dictates that the controls be analyzed before sample results.

#### 8.1.1. Controls

If kit control(s) fail, the test is invalid and needs to be repeated.

Sample data is analyzed and interpreted only after all the kit controls pass.

- NTC (no template control): must NOT have a detectable Ct in the S gene or RNase P reactions. If this control has a detectable Ct in any of the reaction wells, this indicates contamination of the PCR run and it is considered invalid and must be repeated.
- PC**: Ct ≤ 38 for S gene L452R target (FAM™)
- Internal process control (RNase P): all human samples should exhibit Ct values Ct ≤ 38, thus indicating the presence of the human RNase P gene.

Failure to detect RNase P in any specimens may indicate:

- Improper extraction of nucleic acid from clinical materials resulting in loss of RNA and/or RNA degradation.
- (RT-) PCR inhibition.
- Absence of sufficient human cellular material due to poor collection or loss of specimen integrity.
- Improper assay set up and execution.
- Reagent or equipment malfunction.

If the RNase P assay does not produce a positive result for specimens, interpret as follows:

- If S gene marker is positive even in the absence of a positive RNase P, the result should be considered valid. It is possible that some samples may fail to exhibit RNase P amplification curves due to low copy numbers in the original sample. A negative RNase P signal does not preclude the presence of SARS-CoV-2 RNA in a specimen.
- If the S gene marker AND RNase P are negative for the specimen, the result should be considered invalid for the specimen. If residual specimen is available, repeat the extraction procedure and repeat the test. If all markers remain negative after re-test, report the results as invalid and a new specimen should be collected if possible.

#### 8.1.2. Samples

If test run is valid interpretation of sample results is as follows:

- Positive (POS): Ct ≤ 38
- Negative (neg): Not detected (ND) or Ct > 38

The probe for detection of the S gene L452R mutation is labelled with FAM™. The probe for detection of the RNase P gene is labelled with CY5™. Please refer to the following table for interpretation of different signal combinations.

| Detection Channel |                      | Interpretation of Results   |
|-------------------|----------------------|---|
| FAM™ (L452R)      | CY5™ (RNase P)       |   |
| Positive          | Positive or Negative | S gene L452R mutation detected  |
| Negative          | Positive             | S gene L452R mutation not detected.<br>This result is an indication of the presence of the L452 wildtype sequence due to the positive result in the preceding SARS-CoV-2 RT-PCR screening assay |
| Negative          | Negative             | PCR inhibition or reagent failure.<br>A diagnostic statement must not be made.<br>The RT-PCR run should be repeated or a new sample must be analyzed.   |

## 9. LIMITATIONS OF THE PROCEDURE

Very low amounts of starting material, e.g. weak positive samples (high Ct values in the previous diagnostic SARS-CoV-2 RT-PCR), can result in GSD NovaType Select L452R SARS-CoV-2 test runs with negative signals. However, such results do not exclude the presence of SARS-CoV-2 RNA and therefore do not exclude the presence of a L452R mutation.

## 10. TRADEMARKS AND DISCLAIMERS

---

Registered names, trademarks, etc. used in this document are to be considered protected by law even if not specifically marked as such.

## 11. PRECAUTIONS AND WARNINGS

---

- **FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.**
- All human samples should be regarded and handled as potentially infectious.
- Do not interchange reagents of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Wear disposable powder-free gloves, a laboratory coat and eye protection when handling specimens.
- Always use DNase/RNase-free disposable reaction tubes and pipette tips with aerosol barriers.
- Avoid microbial and nuclease (DNase/RNase) contamination of the specimen and the components of the kit.
- In order to avoid contamination of working space with nucleic acids, reaction tubes/plates should not be opened after amplification.
- RT-PCR is highly sensitive to nucleic acid contamination. Therefore, positive / potentially positive material must be stored separately from all other components of the kit.
- Dedicate supplies and equipment to the separate working areas and do not move them from one area to another.
- This assay must not be used on the specimen directly.
- Prior to using this assay the nucleic acid has to be extracted with suitable extraction methods from the original specimen.
- Since ethanol is a strong Real-Time PCR inhibitor, it is necessary to completely eliminate it prior to the elution of the nucleic acid during extraction. If using spin columns with washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step of 10 min at approximately 17000 x g (~ 13000 rpm) before eluting the RNA. For this additional centrifugation step, use a new collection tube.
- The result of this RT-PCR kit may be influenced by potential mutations in the genome of the pathogen if they are located in the primer / probe binding region. Underestimation and/or failure to detect the pathogen may occur.
- RT-PCR inhibitors may also elicit underestimation, false negative results or invalid runs. Therefore, only use nucleic acids extraction kits, which remove RT-PCR inhibitors and which are dedicated for downstream RT-PCR processes.
- The Real-Time PCR is only designed for qualified personnel who are familiar with good laboratory practice and trained in RT-PCR.
- For further internal quality control each laboratory should additionally use known samples.

### 11.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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European Centre for Disease Prevention and Control. SARS-CoV-2 - increased circulation of variants of concern and vaccine rollout in the EU/EEA, 14th update – 15 February 2021. ECDC: Stockholm; 2021.

Robert Koch-Institut (2021): RKI - Coronavirus SARS-CoV-2 - SARS-CoV-2: Virologische Basisdaten sowie Virusvarianten. Available online at [https://www.rki.de/DE/Content/InfAZ/N/Neuartiges\\_Coronavirus/Virologische\\_Basisdaten.html](https://www.rki.de/DE/Content/InfAZ/N/Neuartiges_Coronavirus/Virologische_Basisdaten.html), updated on 2/23/2021, checked on 3/16/2021.

Volz, Erik; Mishra, Swapnil; Chand, Meera; Barrett, Jeffrey C.; Johnson, Robert; Geidelberg, Lily et al. (2021): Transmission of SARS-CoV-2 Lineage B.1.1.7 in England: Insights from linking epidemiological and genetic data. In *medRxiv*, 2020.12.30.20249034. DOI: 10.1101/2020.12.30.20249034.

## **ABBREVIATIONS / ABKÜRZUNGEN / ABRÉVIATIONS / ABBREVIAZIONI / ABREVIACIONES / ABREVIATURAS**

|       |  |
|-------|--|
| Ct/Cq | Threshold cycle (Ct) or quantification cycle (Cq) is the cycle number at which the sample's reaction curve intersects the threshold line |
| NAT   | Nucleic acid Amplification Techniques  |
| NTC   | No Template Control  |

**SYMBOLS KEY / SYMBOLSCHLÜSSEL / EXPLICATION DES SYMBOLES / LEGENDA / SIMBOLOS / TABELA DE SIMBOLOS**

|   |   |
|---|---|
|    | Manufactured by / Hergestellt von / Fabriqué par / Prodotto da / Fabricado por / Fabricado por  |
| <b>RUO</b>  | FOR RESEARCH USE ONLY. Not for use in diagnostic procedures. / NUR FÜR FORSCHUNGSZWECKE. Nicht zur Verwendung in diagnostischen Verfahren. / UNIQUEMENT POUR LA RECHERCHE. Ne pas utiliser dans les procédures de diagnostic. / SOLO DI RICERCA. Non per l'uso in procedure diagnostiche. / PARA USO EXCLUSIVO EN INVESTIGACIÓN. No debe utilizarse en procedimientos de diagnóstico. / ESCLUSIVAMENTE PARA A INVESTIGAÇÃO. Não para uso em procedimentos de diagnóstico. |
| <b>LOT</b>  | Lot Number / Chargenbezeichnung / Numéro de lot / Lotto / Número de lote / Número de lote   |
|    | Expiration Date / Verfallsdatum / Date de péremption / Scadenza / Fecha de caducidad / Data de Validade   |
|    | Storage Temperature / Lagertemperatur / Température de conservation / Temperatura di conservazione / Temperatura de almacenamiento / Temperatura de Armazenamento   |
|    | Protect from Light / Vor Licht schützen / Protéger de la Lumière / Proteggere dalla Luce / Proteger de la Luz / Proteger da Luz   |
| <b>REF</b>  | Product Number / Produktnummer / Numéro de produit / Numero del prodotto / Número del producto / Número do produto  |
|    | Consult Instructions for Use / Arbeitsanleitung beachten / Consulter la notice d'utilisation / Consultare le istruzioni per l'uso / Consulte las Instrucciones de Uso / Consultar as Instruções de Utilização   |
| <b>E-MIX</b>  | RT-PCR Enzyme Mix / RT-PCR Enzym-Mix / RT-PCR Mix d'Enzymes / RT-PCR Mix di Enzima / Mescla de Enzimas RT-PCR / Mistura Enzimática RT-PCR   |
| <b>PP L452R</b>   | Primer-Probe-Mix<br>(for detection of mutation L452R and for detection of RNase P gene / zum Nachweis der Mutation L452 und des RNase P Gens / pour la detection de la mutation L452R et pour la détection de la RNase P / per il rilevamento della mutazione L452R e per il rilevamento del gene RNase P / para la detección de la mutación L452R y para la detección del gen RNasa P / para detecção da mutação L452R e para detecção do gene RNase P)                  |
| <b>NFW</b>  | Nuclease free water / Nukleasefreies Wasser / Eau sans nuclease / Acqua priva di nucleosi / Agua libre de nucleasas / Água isenta de Nuclease   |
| <b>PC</b>   | Positive Control / Positivkontrolle / Contrôle positif / Controllo positivo / Control positivo / Controle positivo  |
|  | Contains sufficient for "n" tests / Ausreichend für "n" Tests / Contenu suffisant pour "n" tests / Contenuto sufficiente per "n" saggi / Contenido suficiente para "n" tests / Conteúdo suficiente para "n" testes  |

**NovaTec Immundiagnostica GmbH**

Waldstraße 23 A6  
63128 Dietzenbach, Germany

Tel.: +49 (0) 6074-48760  
Email: info@NovaTec-ID.com  
Internet: www.NovaTec-ID.com

Fax: +49 (0) 6074-487629

PCOVL452Rx\_RUO\_rev01\_20210518