GSD NovaType II SARS-CoV-2

RT-PCR

RUO

FOR RESEARCH USE ONLY
Not for use in diagnostic procedures

English	2
Deutsch	8
Bibliography / Literatur / Bibliographie / Bibliografia / Bibliografía / Bibliografia	15
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviaturas	15
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos	16

Product Number:	PCOV6083T	(1 x 96 Determinations)
	PCOV6084T	(4 x 96 Determinations)

ENGLISH

1. INTRODUCTION

End of 2019, a novel respiratory disease emerged in the city of Wuhan, Hubei Province of the People's Republic of China, and soon spread rapidly within the country and worldwide. The causative agent was identified as severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). SARS-CoV-2 (2019-nCoV), like the closely related SARS coronavirus (SARS-CoV), belongs to the genus Betacoronavirus within the family of coronaviruses. The zoonotic reservoir of the virus appears to be bats.

Coronaviruses are enveloped, positive single-stranded large RNA viruses that infect humans, but also a wide range of animals. The common human coronaviruses NL63, 229E, OC43 and HKU1 are widespread especially throughout the winter months. They are responsible for up to one third of all acute respiratory diseases, typically with mild symptoms (common cold). More than 80 % of the adult population have antibodies against human coronaviruses. The immunity from previous infections lasts only for a short period of time. Therefore, reinfections with the same pathogen are possible just after one year.

SARS-CoV-2 is predominantly transmitted by droplet infection via coughing or sneezing and through close contact with infected patients. In theory, smear infection and infection through the conjunctiva of the eyes are also possible. The incubation period is in the median 5–6 days (and up to 14 days maximum).

The clinical manifestations of SARS-CoV-2-related COVID-19 disease include fever, cough, respiratory problems and fatigue. In most patients the infection manifests with symptoms of a mild febrile illness with irregular lung infiltrates.

The initial clinical sign of COVID-19 which allowed case detection was pneumonia. But it turned out that the course of the disease is non-specific and varies widely, from asymptomatic courses to severe pneumonia with lung failure and death. However, based on current knowledge, around 80 % of the illnesses are mild to moderate.

Although severe courses of the disease also occur in younger patients and people without previous illness, the following groups of people have an increased risk of serious forms of the disease: elderly people (with a steadily increasing risk from around 50-60 years of age), smokers and people with certain diseases of the cardiovascular system or the lungs, patients with chronic liver diseases, diabetes mellitus, cancer, or patients with a weakened immune system (e.g. due to immune deficiencies or by taking drugs that suppress the immune system).

Treatment focuses on supportive measures depending on the severity of the clinical picture (e.g. oxygen administration, fluid balance management, etc.) as well as the treatment of relevant underlying diseases. Various specific therapeutic approaches (directly antiviral effective, immunomodulatory effective) have been and are being investigated in studies during the course of the SARS-CoV-2 pandemic.

The SARS-CoV-2 trimeric spike glycoprotein (S) is exposed on the surface of virions and mediates viral entry into host cells. The S protein constitutes the primary target of all current leading vaccine candidates. This large, exposed protein is readily targeted by neutralizing antibodies, which indirectly creates selective pressure for the emergence of evasion mutations. The propensity for S to mutate may limit its future use in serological assays and vaccines, as antibodies directed against the current variant may not bind emerging mutated epitopes.

While many countries are currently experiencing a decline in overall infections, most probably due to the impact of tightened non-pharmaceutical interventions, the introduction and increased spread of new SARS-CoV-2 variants first identified in the United Kingdom (B.1.1.7), South Africa (B.1.351), and Brazil (P.1) has raised concerns.

B.1.1.7, also known as VOC 202012/01 (Nextstrain clade 20B, GISAID clade GR and PANGO lineage B.1.1.7), is defined by multiple spike protein changes (deletion 69-70, deletion 144, amino acid change N501Y, A570D, D614G, P681H, T716I, S982A, D1118H) as well as by mutations in other genomic regions. In addition, genomically confirmed B.1.1.7 cases with an additional mutation (E484K) have been reported. This mutation is also carried by the B.1.351 and P.1 variants.

Among others, epidemiological and contact tracing data suggest an increased reproductive number as well as easier transmissibility for the B.1.1.7 variant. This may be due to higher receptor affinity through the N501Y polymorphism. According to current knowledge, the spike protein is the most important target for the action of neutralizing antibodies. Polymorphisms in this protein could therefore affect the effectiveness of the response to vaccination. However, preliminary *in vitro* studies indicate that the licensed mRNA vaccines are also effective against lineage B.1.1.7 viruses.

The B.1.351 (20H/501Y.V2) variant, first identified in South Africa (Nextstrain clade 20C, GISAID clade GH, and PANGO lineage B.1.351) is defined by multiple spike protein changes present in all viruses in the cluster (D80A, D215G, E484K, N501Y and A701V), and more recently collected viruses have additional changes. Three of the changes (K417N, E484K, and N501Y) are located within the receptor-binding domain (RBD).

Increased transmissibility is also discussed for the B.1.351 variant; evidence suggests an increase in ACE2 (Angiotensin-converting enzyme 2) receptor affinity when the E484K and N501Y polymorphisms are combined. The K417N and E484K polymorphisms decrease sensitivity to neutralizing antibodies, which may indicate reduced efficacy of the immune response after infection and vaccination, respectively.

Variant P.1 (20J/501Y.V3, Nextstrain clade 20B, GISAID clade GR and PANGO lineage P.1), first reported by Japan in returning travelers from Brazil, and then later in Brazil has 11 amino acid changes in the spike protein compared to its ancestral lineage B.1.1.28, three of which are located in the receptor-binding domain. The full set of spike protein changes for the variant are L18F, T20N, P26S, D138Y, R190S, K417T, E484K, N501Y, H655Y, T1027I, and V1176F. The P.1 variant resembles the variant described from South Africa in certain RBD key positions (K417, E484, N501). Therefore, increased transmissibility and reduced efficacy of neutralizing antibodies are also discussed for this variant.

Species	Disease	Symptoms e.g.	Transmission route
SARS-CoV-2 (severe acute respiratory syndrome coronavirus 2)	COVID-19	the course of the disease is unspecific, diverse and varies greatly, from asymptomatic courses to severe pneumonia with lung failure and death	primary mode of transmission: droplet infection; smear infections and infections via the conjunctiva of the eyes are theoretically possible

Infection or presence of pathogen may be identified by

- Nucleic acid Amplification Techniques (NAT): e.g. RT-PCR
- Serology: detection of antibodies by e.g. ELISA

2. INTENDED USE (FOR RESEARCH USE ONLY!)

The GSD NovaType II SARS-CoV-2 is a research use only (RUO) real-time reverse transcription PCR (RT-PCR) assay designed for the simultaneous qualitative detection and discrimination of the SARS-CoV-2 (Severe acute respiratory syndrome coronavirus 2) S gene variants wildtype, and mutations K417N, E484K and N501Y. Sample material should be genomic RNA extracted from human upper respiratory (nasal wash/smear, nasopharyngeal wash/smear, oropharyngeal swab) specimen types that have previously tested positive for SARS-CoV-2 RNA by diagnostic RT-PCR methods e.g. GSD NovaPrime® SARS-CoV-2 (COVID-19). The GSD NovaType II SARS-CoV-2 is expressly not intended for diagnosis of SARS-CoV-2 infection and not for use in diagnostic procedures. It is only designed to aid investigations related to the prevalence and frequency of virus variants.

3. PRINCIPLE OF THE ASSAY

The qualitative determination of specific RNA is based on Real-Time reverse-transcription Polymerase Chain Reaction (RT-PCR) technology. The kit contains specific primers and probes labelled with fluorescent reporter and quencher dyes for amplification and simultaneous detection of specific RNA sequences which represent **specific SARS-CoV-2 S gene variants**.

The gene of interest specific probes are labelled with the fluorophores HEX™, FAM™ and ATTO647N (alternative to CY5™) thereby allowing parallel discrimination between S gene wildtype and mutations N501Y, E484K and K417N. The internal process control Ribonuclease P (RNase P) is labelled with the fluorophore ATTO425. The human RNase P gene primer and probe set serves as an internal positive control to monitor sample quality, RNA extraction, and for detection of inhibitors of the PCR reaction.

The GSD NovaType II SARS-CoV-2 RT-PCR has been validated for following real-time PCR platforms:

- AriaMx™ (Agilent Technologies)
- AriaDx™ (Agilent Technologies)

4. MATERIALS

4.1. Reagents Provided

Cap	Symbol	Component	PCOV6083T (PCOV6084T)	
			Volume per Vial [µL]	Number of Vials
green	E-MIX	RT-PCR Enzyme Mix (reverse transcriptase, hot-start DNA polymerase, RNase inhibitor, nucleotides, magnesium, stabilizers and buffer)	550	1 (4)
blue	PP VAR	Primer-Probe-Mix (for discrimination of variants)	215	1 (4)
yellow	PP IPC	Primer-Probe-Mix (for detection of RNaseP)	110	1 (4)
red	PC	Positive Control (plasmid DNA representing SARS-CoV-2 S gene variants)	150	1 (4)
transparent	NFW	Nuclease Free Water	500	1 (4)

4.2. Materials and Equipment needed, but not provided

- Deionized nuclease-free and nucleic-acid free water
- Real-Time PCR instrument, equipped with suitable set of optical cartridges (for already validated instruments refer to 3. PRINCIPLE OF THE ASSAY; for dyes used refer to 7.2.)
Alternative Real-Time PCR instruments might also be appropriate. Their suitability for use with the GSD NovaType II SARS-CoV-2 has to be validated by the user.
- Appropriate Real-Time PCR consumables (e.g. disposable tubes, reaction plates, corresponding optical closing materials)
- Benchtop microcentrifuge
- Centrifuge with a rotor for microtiter plates
- Vortex mixer
- Adjustable pipettes in relation to reaction setup
- Disposable DNase/RNase free pipette tips with filters
- Disposable powder-free gloves

5. STABILITY AND STORAGE

The GSD NovaType II SARS-CoV-2 kit is shipped on dry ice and all components should arrive frozen.

- All components have to be stored between -30 °C and -15 °C immediately after arrival.
- Repeated freeze thaw cycles (more than two) of reagents should be avoided, since this might affect the performance of the kit. Reagents should be frozen in aliquots if they are used intermittently.
- Keep unfrozen storage (e.g. storage on ice) as short as possible.
- Keep the **E-MIX**, the **PP VAR** and the **PP IPC** in the freezer, until you are ready to use it.
- Protect the **E-MIX**, the **PP VAR** and the **PP IPC** from light.

6. SAMPLE PREPARATION

- Extracted RNA or total nucleic acid extracted from human respiratory specimen types (nasal lavage/smear, nasopharyngeal lavage/smear, oropharyngeal swab, bronchoalveolar lavage) previously tested positive for SARS-CoV-2 RNA by RT-PCR methods e.g. GSD NovaPrime® SARS-CoV-2 (COVID-19) is the starting material for the GSD NovaType II SARS-CoV-2.
- The quality of the extracted RNA has a crucial effect on the performance of the entire RT-PCR test system. Make sure that the nucleic acid extraction method is compatible with Real-Time PCR technology.

7. ASSAY PROCEDURE

7.1. Reaction Setup

- Please read the instructions for use carefully **before** performing the assay. Reliability of results depends on following strictly the instructions for use.
- Before use make sure that all samples and reagents are thawed completely, mixed by up and down pipetting or vortexing and centrifuged briefly.
- Important: **E-MIX** is viscous. Briefly spin down the tube to ensure material has not lodged in the cap or side of tube. Be sure to pipette and dispense carefully, and use pipette tips suitable for pipetting viscous liquids.
- The use of **NFW** as no template control (NTC) is highly recommended.
- Define the positions (wells) for samples and controls (**PC** or **NFW**) on the plate.
- Depending on the RT-PCR instrument used and the available detection channels, RNase P can be used as an internal process control.

Important: If the Real-Time PCR instrument is not equipped with an ATTO425 channel, RNase P amplification must not be performed, i.e. the reaction setup must be pipetted **without** the **PP IPC** Primer-Probe-Mix, otherwise unspecific signals may occur in the other detection channels.

Reaction Setup		
Component	Including amplification of RNase P (ATTO425 channel required)	without amplification of RNase P
E-MIX	5 µL	5 µL
PP VAR	2 µL	2 µL
PP IPC	1 µL	-
NFW	-	1 µL
Sample or PC or NFW	12 µL	12 µL
Total Volume	20 µL	20 µL

- Close the optical reaction plate with corresponding optical closing material.
- Centrifuge the optical reaction plate in a centrifuge with a rotor for microtiter plates for 60 seconds at approximately 1000 x g (~ 3000 rpm).

7.2. Programming the Real-Time PCR Instrument

Regarding setup and programming of the Real-Time PCR instrument, please use the manual of the respective instrument.

For detailed programming instructions regarding the use of the GSD NovaType II SARS-CoV-2 RT-PCR on specific Real-Time PCR instruments please contact NovaTec Immundiagnostica GmbH.

RT-PCR Run Settings	
Reaction Volume	20 µL
Passive Reference	None

Fluorescent Detectors / Fluorophores

Detection of the amplified viral nucleic acid fragment is performed in the detection channels FAM™ (S gene mutation N501Y), HEX™ (S gene N501 wildtype and mutation K417N) and CY5™ (S gene mutation E484K; fluorophore: ATTO647) with MGBEQ Quencher.

Detection of the amplified internal process control RNase P is performed in the detection channel ATTO425 with BHQ1 Quencher.

Detection	Wildtype	B.1.1.7 (501Y.V1)	B.1.351 (501Y.V2)	P.1 (501Y.V.3)
K417N mutation	-	-	HEX™	-
E484K mutation	-	-	CY5™	CY5™
N501Y mutation	HEX™ (wildtype)	FAM™	FAM™	FAM™
RNaseP*	ATTO425	ATTO425	ATTO425	ATTO425

* Detection of RNase P only possible if the cycler used is equipped with the appropriate detection channel (ATTO425) and the **PP|IPC** Primer-Probe-Mix is added to the reaction setup.

Temperature Profile and Data Collection

No. of Cycles	Temperature	Time (min)	Data Collection
1	50 °C	10:00	-
1	95 °C	03:00	-
40	95 °C	00:10	-
	60 °C	00:30	Fluorescence measurement at the end of every cycle

Before starting the test run, please check the settings for cycles, temperature and time.

8. RESULTS

Data analysis should be performed with the software of the used real-time PCR device according to manufacturer's instructions.

Analysis settings:

Setting	Recommendation
Threshold	Enter a value for the threshold so that the threshold is: <ul style="list-style-type: none"> Above the background Below the plateau and linear regions of the amplification curve Within the exponential phase of the amplification curve
Baseline	Select start and end cycle values so that the baseline ends before significant fluorescence is detected.

8.1. Interpretation of Results

Test run is valid only if RT-PCR run complete.

The GSD NovaType II SARS-CoV-2 test protocol dictates that the controls be analyzed before sample results.

8.1.1. Controls

If kit control(s) fail, the test is invalid and needs to be repeated.

Sample data is analyzed and interpreted only after all the kit controls pass.

- NTC (no template control): must NOT have a detectable Ct in the S gene or RNase P reactions. If this control has a detectable Ct in any of the reaction wells, this indicates contamination of the PCR run and it is considered invalid and must be repeated.
- PC**: Ct ≤ 38 for all S gene targets (HEX™, FAM™ and Cy5™)
- Internal process control (RNase P): all human samples should exhibit Ct values Ct ≤ 38, thus indicating the presence of the human RNase P gene.

Failure to detect RNase P in any specimens may indicate:

- Improper extraction of nucleic acid from clinical materials resulting in loss of RNA and/or RNA degradation.
- (RT-) PCR inhibition.
- Absence of sufficient human cellular material due to poor collection or loss of specimen integrity.
- Improper assay set up and execution.
- Reagent or equipment malfunction.

If the RNase P assay does not produce a positive result for specimens, interpret as follows:

- If S gene markers are positive even in the absence of a positive RNase P, the result should be considered valid. It is possible that some samples may fail to exhibit RNase P growth curves due to low cell numbers in the original sample. A negative RNase P signal does not preclude the presence of SARS-CoV-2 RNA in a specimen.
- If all S gene markers AND RNase P are negative for the specimen, the result should be considered invalid for the specimen. If residual specimen is available, repeat the extraction procedure and repeat the test. If all markers remain negative after re-test, report the results as invalid and a new specimen should be collected if possible.

8.1.2. Samples

If test run is valid interpretation of sample results is as follows:

- Positive (POS): Ct ≤ 38
- Negative (neg): Not detected (ND) or Ct > 38

SARS-CoV-2		Detection Channel			
S gene	Variant	HEX™	FAM™	CY5™	ATTO425* (internal process control)
N501	wildtype	POS	neg	neg	POS
N501Y	B.1.1.7 (501Y.V1)	neg	POS	neg	POS
K417N, E484K, N501Y	B.1.351 (501Y.V2)	POS	POS	POS	POS
E484K, N501Y	P.1 (501Y.V.3)	neg	POS	POS	POS

* The RNase P gene is used as an internal process control because many copies of it exist in the human genome, and it is readily detectable. RNase P is derived from the human cells that are present in every sample used. It is assumed that if human nucleic acids were extracted to detect the human gene RNase P, viral nucleic acids were also successfully extracted. The RNase P target is also amplified as a quality control for the extraction method and to corroborate the absence of PCR-inhibitors in the sample.

Discrimination of virus variants with respect to S gene codons 417, 484 and 501 is only possible for samples

- that have already been pretested positive for the presence of SARS-CoV-2 RNA and
- whose Ct value(s) in the corresponding S gene channel(s) is equal or below 38 (Ct ≤ 38)

Note on variant B.1.1.7:

1. Other spike protein changes, e.g. deletion 69-70 do not affect the present RT-PCR, as these changes are located far from the detected amplicons.
2. With the present setup, B.1.1.7 variants carrying an additional E484K mutation cannot be distinguished from the P.1 variant, as both show signals in the FAM™ and Cy5™ channels.

9. LIMITATIONS OF THE PROCEDURE

Very low amounts of starting material, e.g. weak positive samples (high Ct values in the previous diagnostic SARS-CoV-2 RT-PCR), can result in GSD NovaType II SARS-CoV-2 test runs with negative signals. However, such results do not exclude the presence of SARS-CoV-2 RNA.

10. TRADEMARKS AND DISCLAIMERS

Registered names, trademarks, etc. used in this document are to be considered protected by law even if not specifically marked as such.

11. PRECAUTIONS AND WARNINGS

- **FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.**
- All human samples should be regarded and handled as potentially infectious.
- Do not interchange reagents of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Wear disposable powder-free gloves, a laboratory coat and eye protection when handling specimens.
- Always use DNase/RNase-free disposable reaction tubes and pipette tips with aerosol barriers.
- Avoid microbial and nuclease (DNase/RNase) contamination of the specimen and the components of the kit.
- In order to avoid contamination of working space with nucleic acids, reaction tubes/plates should not be opened after amplification.
- RT-PCR is highly sensitive to nucleic acid contamination. Therefore, positive / potentially positive material must be stored separately from all other components of the kit.
- Dedicate supplies and equipment to the separate working areas and do not move them from one area to another.
- This assay must not be used on the specimen directly.
- Prior to using this assay the nucleic acid has to be extracted with suitable extraction methods from the original specimen.
- Since ethanol is a strong Real-Time PCR inhibitor, it is necessary to completely eliminate it prior to the elution of the nucleic acid during extraction. If using spin columns with washing buffers containing ethanol, it is highly recommended to perform an additional centrifugation step of 10 min at approximately 17000 x g (~ 13000 rpm) before eluting the RNA. For this additional centrifugation step, use a new collection tube.
- The result of this RT-PCR kit may be influenced by potential mutations in the genome of the pathogen if they are located in the primer / probe binding region. Underestimation and/or failure to detect the pathogen may occur.
- RT-PCR inhibitors may also elicit underestimation, false negative results or invalid runs. Therefore, only use nucleic acids extraction kits, which remove RT-PCR inhibitors and which are dedicated for downstream RT-PCR processes.
- The Real-Time PCR is only designed for qualified personnel who are familiar with good laboratory practice and trained in RT-PCR.
- For further internal quality control each laboratory should additionally use known samples.

11.1. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

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