


EUROFINS GENOMICS

The DNA Universe



GRAS-Di[®] for Marker discovery & Genotyping

- **GRAS-Di[®]** (**G**enotyping by **R**andom **A**mplicon **S**equencing-**D**irect) was developed by **TOYOTA Motor Corporation**
- GRAS-Di[®] amplifies the random regions throughout the genome for NextGeneration Sequencing.
- Skim based sequencing of whole genome covered by the amplicons through 2 step PCR with random primers, Non-targeted PCR-based GBS.
- Sequencing data is analyzed for presence / absence of amplicon markers and optionally SNPs calling.




GRAS-Di[®] vs Whole Genome Sequence(WGS)

- **GRAS-Di[®]**: Sequence reads are uniformly mapped to the genome **at intervals**.
- **WGS**: Sequence reads are uniformly mapped to the genome **without gaps**.




Increasing of sequence data per samples.




Amplicons covering on the chromosomes

- Mapping of Amplicons onto the reference sequence shows the uniform distribution of Markers over the genome



GRAS amplicons were distributed over genome, uniformly.

- More cost effective than Whole Genome deep Sequencing



- Economical & Ecological, simple method that tolerates low quality DNA or slightly fragmented DNA (>100ng)
- Applicable to both natural and segregated populations, all without the need for a reference genome
- Obtains fewer missing data than RAD-seq

Applicable for Genetic diversity studies, Polymorphism detections (SNPs or InDels or Haplotypes).
Increased usage in the field of Phylogenetic and population structure studies.

SNP Discovery

- Develop SNP markers for translation to a preferred genotyping system
- High quality alternative to GBS/WGS, RNA-Seq, RAD-Seq, ddRAD-Seq (GBS)

High Throughput Genotyping

- Discover and genotype massive amounts of genetic variants for Marker Selection
- High-density Linkage Mapping & QTL Mapping, Chromosome Mapping, Phylogenetic & Population structure studies

Others

- Potential of large number of polymorphisms used as molecular markers for genetic analysis
- Possible to obtain genome-wide co-dominant markers Identify many markers covering all chromosomes even in a population with small genetic variation

GRAS-Di[®] has been implemented successfully in **over 120 different species** including crops, forestry, livestock and aquaculture.

The technology is also applicable to highly polyploid species.
Please refer the Bibliography for further details.



Agriculture/Plant



Forestry/Fruit tree



Livestock/Animal



Fishery/Aquaculture

Rice, Wheat, Corn, Sugar cane, Soy bean, Groundnut, Green pea, Tomato, Eggplant, Green pepper, Paprika, Potato, Sweet potato, Cabbage, Lettuce, Chinese cabbage, Japanese white radish, Broccoli, Cauliflower, Green onions, Onions, Cucumber, Pumpkin, Spinach, Carrot, Burdock, Strawberry, Melon, Watermelon, Sub-clover, Horse gram, *Lotus japonicus*, Apple, Peach, Pear, Orange, Persimmon, Grapes, Fig, Kiwi, *Citrus junos*, *Citrus sudachi*, Prune, *Chrysanthemum*, Cedar, Cherry tree, Human, Cattle, Pigs, Chicken, Horse, Rabbit, Hamster, Guinea pigs, Rat, Tuna, Shiitake mushroom, Hen of the Woods

- Raw sequence data in FASTQ format
- Genotyping data by TOYOTA software with Dominant markers

Analysis	Methodology	Deliverables	Requirements
Toyota software	Genotyping based on presence or absence of amplicon reads.	GRAS-Di_RESULT.csv	Information about population (natural population or segregating population (including parent and progeny samples))
		GRAS-Di_RESULT_codominant.csv	Segregating population, segregation ratio of the genotypes in the segregating population (eg. F2 population 1:2:1), diploid, homozygous parents
		GRAS-Di_RESULT_mapping.csv	reference genome
SNP calling	Mapping and calling SNPs	*.bam	reference genome
		*.vcf	reference genome

(Optionally)

- Genotyping result with Co-dominant markers
 - * Available only in case of diploids, segregated population & homozygous parents, more than 40 samples of F2 generation.
- Mapping result of dominant markers to the reference genome
 - * Available only for the species whose reference genome is available
- Mapping & SNP call result
 - (1) Mapping result (BAM format)
 - (2) SNP call result (VCF format)


Appendix; Result of TOYOTA software

▽ “GRAS-Di” is provided the co-dominant marker in diploid segregation population.


	A	B	C	seq1	seq2	E			
	PARENT					length			
1									
2	C1	TEST0000082	AAGGGCAACCAGCA	CCACAGATAGCA	ACCAATTAACCA	ACCCGAGCC	CGGTAGACGTGCTGGGCTCGGGTTGGTTAATTGGTTTGCATCTGTGGG	66	SNP
3	C2	TEST0000085	AACGGCAACCAGCA	CCACAGATAGCA	ACCAATTAACCA	ACCCGAGCC	CGGTAGACGTGCTGGGCTCGGGTTGGTTAATTGGTTTGCATCTGTGGG	66	SNP
4	C2	TEST0000132	AAGAGAGCGCAAAAGTCAAAAACCTAGCGTTCCGCT	TTCCACA	ACTGCCA	CGAATCACATCCTCCTCGACTGTGTGTTTGCCTCGACAGGTCTGGCTCCGCGT	999	SNP	
5	C1	TEST0000135	AAGAGAGCGCGGAGTCAAAAACCTAGCGTTCCGCT	TTCCACA	ACTGCCA	CGAATCACATCCTCCTCGACTGTGTGTTTGCCTCGACAGGTCTGGCTCCGCGT	999	SNP	
6	C1	TEST0000139	AAGAGGAGGCCAC	TAATTGGTGTCACTGTGCTGCCCTCATCACTGACGTACA	AGCGGGTCTGCAACAAATAGATCTAGAGGTCTAGGCCCCCTCCCCGGTTGGAT	114	Haplo-type		
7	C2	TEST0000141	AAGAGGAGGCCAC	TAATTGGTGTCACTGTGCTGCCCTCATCACTGACGTACA	AGCGGGTCTGCAACAAATAGATCTAGAGGTCTAGGCCCCCTCCCCGGTTGGAT	114	Haplo-type		
10	C1	TEST0000240	ACAAACCAAGCAAAGCGATATC	AAATCCTAGCCATTAATCGATGTCGGG	CTCAAAAGTAGGAAAAAGGCAAAAGAAATGACCAGAACAGTGGTCCCCGACAT	96	SNP		
11	C2	TEST0000241	ACAAACCAAGCAAAGCGATATC	AAATCCTAGCCATTAATCGATGTCGGG	CTCAAAAGTAGGAAAAAGGCAAAAGAAATGACCAGAACAGTGGTCCCCGACAT	96	SNP		
12	C1	TEST0000246	ACAAACTGCBAACTGCTACAACCTCACGSGCAAAAGAGCAATTGAAGACAGT	AGCAGCA	AGCAAGAAACCGGATCCATCCATCCATCCATCCATAGGTAATG	142	SNP		
13	C2	TEST0000247	ACAAACTGCBAACTGCTACAACCTCACGSGCAAAAGAGCAATTGAAGACAGT	AGCAGCA	AGCAAGAAACCGGATCCATCCATCCATCCATCCATAGGTAATG	142	SNP		
14	C2	TEST0000278	ACAACTACAAAGCCAAACATACGACAAGACAAGACGTCCGATATCGACTTCGG	TAATGTCTCTCCTATCAGTATCCGGAGACACCGTAGGGGACTAGCCGTGC	CTATCTCCGAAGTCCGATATCCGACGCTTGTCTTGTCTGATGTTGGCT	110	SNP		
15	C1	TEST0000279	AGATAGGCACGGCTAGTCCCCTACGGTGTCTCCGGATACTGATAGGAG	TAATGTCTCTCCTATCAGTATCCGGAGACACCGTAGGGGACTAGCCGTGC	CTATCTCCGAAGTCCGATATCCGACGCTTGTCTTGTCTGATGTTGGCT	110	SNP		

“GRAS-Di” could be identified markers by all sequence information.

- Genotypes of each locus are determined by GRAS-Di software based on the presence or absence of amplicons (amplicon is considered as locus).
- When amplicon sequences are different between parents, and a reference genome is available, then the amplicons are available for SNP calling.



Appendix; Amplicon Mapping on Chromosome



id	Seq1	Seq2	ref ID	pos	num mistakes	ref ID	pos	num mistakes	ref ID	Pos	num mistakes	ref ID	pos	num mistakes
Amplicon 1	ATCGTGTTCGGACGGGAATAACTGTGTAATG AAOOGCAATAATGCGTGTACATTTCGGTCGAGT DCAGGAATGGGGGGCGGACTGGGGGACAAAG	TACCGGTTTTGTTGTCGGGGCAGTGGGGGGCCGA TTGCTGGACTCGAOGGAAATGTAAGAGGATATT GGCGTTGATTGACAGTATTATTTCGGTCGGG	Chr 01	43,216,785	4	Chr 07	334,562	4	Chr 12	16,339,644	4	Chr 06	19,715,592	5
Amplicon 2	ATCGGTAGTTAGCTTGGGCTTTGTCGGAGGTAC AAGGGGGACTATAATAATGCTTTTGGGCTAAAGG GGATATAGATGGCTGTACGGTCATCATOTGC	TGAATCTATGGACCTTTGTTGTCGACCTTCGGGGGCG TTGCTAAGTGTACGTTTGTAGGGGCTATAGGGGGCGC AGGTGTTGCGATGTGTACACTGAAAGGGGG	Chr 02	1,155,876	1	Chr 12	25,516,783	2	Chr 10	14,005,983	2	Chr 01	26,208,751	2
Amplicon 3	ACAGGGTCCACAAGGTGAGAGGAGTTAGATAC GACCTATATAGTCTGGGGGCTATTTGGCATAATA AAATGGTCAGGACGGGAGGTT	TGCATAGCATAAGTTGGTTGATGGAGGGCAGTCAAAAG TTGCACTGTAGATCCGACGATTTGGAGGGGATTAATA GGGTTAGTCGGTAACCTGGATGATGATGGG	Chr 03	613,758	1	Chr 03	616,789	1	Chr 03	617,946	2	Chr 05	11,003,782	3
Amplicon 4	ACAGGGTCCACAAGGTGAGAGGAGTTAGATAC GACCTATATAGTCTGGGGGCTATTTGGCATAATA AAATGGTCAGGACGGGAGGTT	TGCATAGCATAAGTTGGTTGATGGAGGGCAGTCAAAAG TTGCACTGTAGATCCGACGATTTGGAGGGGATTAATA GGGTTAGTCGGTAACCTGGATGATGATGGG	-	-	-	-	-	-	-	-	-	-	-	-
Amplicon 5	GGCTAGTTGCTGGATTTGCAAGTGGTTCAAGT TTTTTAATAGTGGCTAAGCAGTACTAATATGAG TGTATGAGGTTAAAGGGGGGCTATGGGGG	ACTCTGCTTTAGCTCACAGATTACGATTTGCTGGC ATGGGGGATAGGGGGGGTTTAAAGCTTCATACAC TGATTTAGTACTGGTTAGGGACTATTAATA	Chr 05	9,994,255	0	-	-	-	-	-	-	-	-	-

Multiple candidates of mapped positions for each amplicon are displayed

- ref ID: Chromosome name
- pos: coordinate position on the chromosome
- num mistakes: Number of mismatch bases with reference genomic sequence

Please Contact to nearest Eurofins Genomics

For APAC Region

**Eurofins Genomics
K.K., Japan**

E-mail:
ngs-jp@eurofins.com



For EU Region

**Eurofins Genomics
Germany GmbH**

E-mail:
projectsales-eu@genomics.eurofinseu.com

